

Options

Service pack 2

New Module Slide

Service Port

Completely Removable Cover



Options

www.liconic.com

Europe: Liconic AG, Liechtenstein

Liconic UK Ltd, Manchester

US & International: Liconic Inc, Woburn MA

SP2

Service pack 2

The new service pack 2 includes several innovative features that improve serviceability of the LiCONiC incubator family. These features provide our units to stay connected to the workstation while monitoring all functionalities of the system. This ensures the units never have to be re-aligned during preventative maintenance or service.

New Module Slide:

All our components are on the new module slide, which allows service engineers slide the whole module out for easy access to every component without removing the incubator from the workstation. With all components on this module, inspection and replacement of parts can be done in minutes. This saves valuable down time and realignment times, which occur with other types of incubators.

Additional Comports (UCI?):

Additional comports allow the end-user to run additional software on the host in parallel to the existing driver software. Monitor all aspects of the incubator, without disconnecting from the work station or upper level software. This function allows the end-user to simultaneously monitor PLC, Temp, RH%, and gassing functions of the unit and store as log files, which can be used for set-up or GLP and QA statistics.

Service Port:

Use the new Service Port for system diagnosis and system performance monitoring. This port may also be used for firmware updates, firmware backups or system configuration using LiCONiC's simple-to-use software tools. This service port will also support Liconic powerful software tools for handling set-up. The service port allows engineers to plug a PC to the unit and perform desired functions and monitor all aspects of the incubator, without disconnecting from the work station or upper level software.

Completely Removable Cover:

The advantage to the removable cover is total access to refrigeration compressor if needed. It opens the work environment for field engineers which will save time and money on repairs.

| <i>SP2 option ...</i> | <i>Order Nr.</i> |
|------------------------------|-------------------------|
| STX40 | 9118 12 47 |
| STX220 | 9123 00 22 |