

<b>Title:</b>	<b>Cooling does not work</b>			<b>ID:</b>
				0063
<b>Date in:</b>	<b>Response:</b>	<b>Model:</b>	<b>Author:</b>	
2006-08-18	2006-08-18	STX1000	CMa	

**Q:**

On an STX1000-DFSA (Deep Freezer) we have +38'C.

**A:**

First set of information

- What is the temperature on the lower display?
- Does the compressor run?
- Check Overpressure Switch
- Does blower inside chamber run
- Is Relay 515 activated
- Is Relay 513 activated

