

Title:	Plate Capacity STX	ID:	
		0167	
Date in:	Response:	Model:	Author:
2003-01-03	2003-01-03	STX40 / STX200	CMa

Q:

I would like some definitive answers as to the capacities of your STX200, STX40 incubators and your plate hotel. In particular, how many 384-well plates can these devices accomodate. Nobody here seems to be quite sure.

A:

To your question: Assuming a plate + lid height of 17mm there

STX40 44 plates

STX200 198 plates (Models 2nd Q. 2003: 220 plates)

Please let me know if the plates are less high.