

Title:	Stacker for Tip Boxes			ID:
				0208
Date in:	Response:	Model:	Author:	
2004-03-11	2004-03-11	-	C <i>Ma</i>	

Q:

Before we make any mistake I want to know the maximum height in the deelwell cassette.

We are going to use :

DeepWell Racks - height = 44.5 mm

Tip Boxes from - height = 50.5 mm

Is this OK for the standard deepwell cassette ?

A:

Pitch is 50mm so its 45mm. For Tip Boxes we have P72, 7 levels that are slightly wider than standard cassettes.