

Accessories – Transfer Stations

SA-Attachment Kit



The SA-Attachment Kit is required only for LiCONiC's STX44-xxSA and STT1000-Kiwi. Normally, the BT-Attachment Kit can be mounted on all LiCONiC incubators. However, the standard mounting position for the BT-Attachment is not available for LiCONiC's STX44-xxSA and STT1000-Kiwi. The kit can be used by customers who have changed the design of their work-station and need a standard solution for the next design generation.

The components used are the most reliable and of the highest quality available in the market. All manufactured components are electro-coated to protect against corrosion and contamination from external environmental conditions.

Specifications

Available for	STX44-xxSA and STT1000-Kiwi
Design	Fixture for Transferstations