

## Accessories – Transfer Stations

### Slide Station

The **function** of the slide-station was designed for the use of non-articulating robotic arms that need additional clearance from the storage units. With its precise stopping point it ensures these devices need to be taught only once for accurate load and retrieval of the storage plates.

The components used are the most reliable and highest quality products available. All of our manufactured components are electro-coated for the prevention of corrosion and contamination from external conditions.



#### Order information

*For STX, LPX Series*

	<u>Order Nr.</u>
Slide transfer station for extended transfer reach motor up (170 mm)	9118 12 28
Slide transfer station for double extended transfer reach (300 mm)	9111 07 81
Custom size slide transfer station (specify transfer reach 270-1000 mm)	9118 12 30

#### Specifications

<b>Plate Positions</b>	1
<b>Sensor</b>	Plate presence sensor
<b>Layout</b>	Straight on
<b>Movement</b>	Perpendicular to / from Unit
<b>Special</b>	To bridge a gab