


Product introduction release

STR 240

The STR is the most versatile incubator in the market. Fits underneath any work-station or liquid handling system. Reduces your bottlenecks by 50%, when using multiple plate transfer nests. You save 30-60K on accessories because you do not have to buy extended or articulating robotic arms. The periscope delivery shovel presents the plate to the top of the work surface. With work space at a premium, the STR will save valuable space for other components above your system and gives you a more aesthetically pleasing compact appearance.



All this without sacrificing price. 240 plate capacities in such a foot print is cutting edge technology. Visit www.liconic.com for more details.