

<b>Title:</b>	<b>Transfer Station Inputs 100 and 101 on MPH</b>			<b>ID:</b>
				0238
<b>Date in:</b>	<b>Response:</b>	<b>Model:</b>	<b>Author:</b>	
2004-09-11	2004-09-11	MPH	CMa	

## Q:

Cosmas, I cannot read 100 or 101, but I know the inputs are working since with out a jumper I am unable to burn in the unit. Could you please look into these issues and get back to me.

## A:

The MPH firmware goes back to the cytomat times. On the MPH

100 is on 007  
101 is on 006