

Title:	Various Plate Heights and Tip Boxes	ID:	
		0249	
Date in:	Response:	Model:	Author:
2004-11-26	2004-11-26	-	CMa

Q:

I have a questions regarding all the incubators. Can they be configured to hold standard plates and tip boxes in the same incubator? I thought it can hold either 1 or the other and that you can not mix different heights in the same unit. Along the same note, tip boxes are slightly taller then deepwell plates. Can we mix all 3 into the same instrument?

A:

Yes, they can handle various plate types and/or tip boxes. Note that very long tips (e.g. 100mm length) may require a special Gate. Please specify maximum height of goods. Please verify also that your user interface in your software allows to work with various pitches. You will have to allow the user to configure the different cassettes.