

Title:	Cleaning a Shaker	ID:	
		<i>0034</i>	
Date in:	Response:	Model:	Author:
<i>2005-06-02</i>	<i>2005-06-02</i>	<i>STX40</i>	<i>CMa</i>

Q:

Do I have to remove the stackers and the shaker for cleaning the incubation chamber of an STX40 Shaker?

A:

We recommend the EasyClean with the Shaker option. This will allow to easily remove the handling with the stackers from the chamber. For cleaning the waterbath the stackers may stay in the system.