

<b>Title:</b>	<b>CO2 Connection</b>	<b>ID:</b>	
		0038	
<b>Date in:</b>	<b>Response:</b>	<b>Model:</b>	<b>Author:</b>
2004-08-02	2004-08-02	-	CMa

## Q:

I have a user at Novartis with some questions regarding their CO2 connection with the StoreX.

- 1) Is it okay to use copper tubing for the connection?
- 2) What should the inlet pressure of CO2 be regulated to?

## A:

- 1) It sure is: The tubing is 4x6mm.
- 2) We recommend ~2 bar pressure. 1..6 is ok