

Plate-Garage

New Plate Mass Storage
for Controlled and Uncontrolled Environment



Over 1.5 Million 0.3 ml Tube Capacity 2.500 - 20.000 MTP Plate Capacity



The new Liconic STL-Series is a family of automated mass storages for 2'500 to 20'000 Microtiter size plates capacity. The STL systems are available with or without controlled environment.

DF: -20°C

DC: 04..25°C 02..80%RH

HC: 04..25°C max. 90%RH

IC: 33..50°C, max. 95%RH

HR: 04..50°C

Liconic Instruments / www.liconic.com

