


Title:	<b>Motor-Card</b>	ID:	126
		14.03.13	Response:
		-	CMA

## Motor Current SMC



Cannel	Reistors (2)	Identification	P.Nbr	Instrument	Application
1	33k	3302		STX1000	Turn
2	15k	1502			Lift
1	22k	2202			Carrousel Right
2	22k	2202			Carrousel Left

Cannel	Reistors (2)	Identification	P.Nbr	Instrument	Application
1	56k	5602		GFX	Rotation
2	39k	3902			Lift
1	56k	5602			x-Axis
2	39k	3902			--