

Title:	Single Phase Power Cables	ID:	
		0097	
Date in:	Response:	Model:	Author:
2010-06-16	2010-06-16	-	CMa

Single Phase Power Cables

Standard incubators like STX-Series, STR-Series etc. and Open Stores like LPX-Series or LPR Series instruments are shipped with a standard power cable.



Depending on the country where the instrument is shipped to, the plug may vary. Commonly used connector types are US (Type B / NEMA), European (Type F / CEE 7/4) or Swiss (Type J / SEV).

For instruments running at higher power ratings there may be an additional high power cable required.



Connectors type supported are US (Type N6/15 / NEMA) and Swiss (Type 23 / SEV). Countries not supported may require local installation of a suitable connector.